

Recumbent Stepper 7.5S



Make rehabilitation accessible for all with the 7.5S stepper, a semi-recumbent designed for wheelchair-bound patients. This stepper aligns with ADA compliance to facilitate accessibility and inclusivity in medical settings.

Patients with even the most limited abilities will benefit from features such as:

- Direct wheelchair access or 360-degree swiveling seat
- Bilateral movement and 20 levels of resistance
- 10-210 steps/min speed range for a variety of patient's needs
- 1"-12.5" adjustable step range

This innovative stepper offers low-impact, full-body exercise to aid in rehabilitation and improves cardiovascular health. Its adjustable settings ensure a safe and personalized experience for individuals of varying abilities.



Recumbent position reduces strain on joints and muscles



Wireless heart rate monitoring capabilities



Direct wheelchair access

Specifications

Braking Device	Eddy current brake with small flywheel for low starting inertia
Drive Train	Heavy-duty, 8-groove poly V belt with high-end linear bearings and guides
Power	100-240 volts AC (standard power supply)
Workload Range	5 to 1000 Watts
Recommended User Height	4'10" - 6'6"
Max User Weight	440 lb / 200 kg
Dimensions	77.6" x 33.5" x 46.5"
Product Weight	305.8 lb
Compliance	CAN/CAS-C22.2 No. 60601-1:14 ANS/AAMI ES60601-1:2005+A2 (R2012) +A1 IEC 60601-1-2:2014 EN 60601-1-2:2015 IEC 60601-1:2015+A1:2012 EN 60601-1:2006+A1:2013+A12:2014 EN ISO 20957:2013 MDD 93/42/EEC Class Im, NB: 0123
Warranty - Commercial (Dues paying)	3 years parts, 3 years labor
Rehab / Light Commercial / Residential (Non-dues)	5 years parts, 1 year labor