

Korebalance® Premiere-22 (shown) or Professional-19



- 22" medical grade touchscreen computer (KB22) or 19" touchscreen computer (KB19)
- Adjustable Monitor Mount Arm
- Individually adjustable safety rails
- Platform has 360° horizontal & 20° vertical range of movement.
- **(KB22) Includes Cognitive, Parkinson's software suite for cognitive testing/training.**
- (KB19) Optional Cognitive, Parkinson's software suite.
- Includes Korebalance® testing/training software suite & games.
- Exercise band attachment points
- Large, stable platform with low step-up height for safe mount/dismount
- Patient weight capacity: 500 lbs. (227 kg)
- Machine weight: 315 lbs. (143 kg)
- Footprint: 44"x58" (112cm x 148cm)
- Max Height: 80" (203cm)
- Power: 100-240 VAC; 50/60 Hz; 500 Watt
- Software available in English, Spanish, French, German & Mandarin
- Frame Color: White
- Warranty: Frame-15 Years Electronics - 2 Years (KB 22) Electronics - 1 Year (KB 19) Wear Items - 1 Year
- Made in the USA

Korebalance® 19 Standard



- 19" touchscreen computer
- Fixed monitor mount arm
- Fixed safety rails
- Platform has 360° horizontal & 20° vertical range of movement.
- Optional Cognitive, Parkinson's software suite for cognitive testing/training.
- Includes Korebalance® testing/training software suite & games.
- Exercise band attachment points
- Large, stable platform with low step-up height for safe mount/dismount.
- Patient weight capacity: 500 lbs. (227 kg)
- Machine weight: 315 lbs. (143 kg)
- Footprint: 44" x 58" (112cm x 148cm)
- Max Height: 80" (203cm)
- Power: 100-240 VAC; 50/60 Hz; 500 Watt
- Software available in English, Spanish, French, German & Mandarin
- Frame Color: White
- Warranty: Frame-15 Years Electronics - 1 Year Wear Items - 1 Year
- Made in the USA

KoreSway Kompact



- 19" OR 22" touchscreen computer with mount allowing both tilt and height adjustments.
- Portable frame with arm rails
- Force Plate Platform with 23" (58.4 cm) diameter
- Optional Cognitive & Parkinson's software suite for cognitive testing/training.
- Includes Korebalance® testing/training software suite & games.
- Patient weight capacity: 330 lbs. (150 kg)
- Total Weight: 112 lbs. (50 kg)
- Max Height: 73" (185.42 cm)
- Width: 40" (101.6 cm)
- Compact and easy to transport
- Software available in English, Spanish, French, German & Mandarin
- Frame Color: White
- Warranty: Frame - 5 Years Electronics – 1 Year
- Made in the USA

KoreSway Mobile or Korebalance Active



- 17.3" Touch Screen Laptop
- Optional Adjustable Laptop Stand and Carry Case.
- **Choice Between:**
 - **KoreSway Mobile** - Force Plate Platform with 23" (58.4 cm) diameter
 - **Korebalance Active** – Variable Stability Platform with 23" (58.4 cm) diameter
- Optional Cognitive & Parkinson's software suite for cognitive testing/training.
- Includes Korebalance® testing/training software suite & games.
- Patient weight capacity: 330 lbs. (150 kg)
- Total weight: 53 lbs. (24 kg)
- Software available in English, Spanish, French, German & Mandarin
- Warranty: 1 Year
- Made in the USA

KoreSway Mini

Korebalance® Trainer



- Variable Stability Platform has 360° horizontal & 20° vertical range of movement
- Platform Diameter: 23" (58.4cm)
- Patient weight capacity: 500 lbs. (227 kg)
- Machine weight: 22 lbs. (10 kg)
- Optional Hand Pump
- Warranty: 1 year
- Made in the USA

FEATURES:	PREMIERE 22	PROFESSIONAL 19	STANDARD 19	KORESWAY COMPACT
ITEM #	KB22 Premiere	KB19 Professional	KB19 Standard	KB COMPACT
COMPUTERIZED BALANCE ASSESSMENT	•	•	•	•
SCREEN SIZE	22"	19"	19"	19" OR 22"
TOUCHSCREEN	•	•	•	•
MEDICAL GRADE TOUCHSCREEN	•	○	○	○
FULLY ADJUSTABLE SCREEN	•	•	○	○
SCREEN ADJUSTABILITY	VERTICAL ADJUSTMENT/TILT	VERTICAL ADJUSTMENT/TILT	TILT	ADJUSTABLE STAND
PLATFORM SIZE	23" DIAMETER	23" DIAMETER	23" DIAMETER	23" DIAMETER
STEP UP HEIGHT	7"	7"	7"	2 ¾ "
PATIENT HAND RAILS	•	•	•	•
INDEPENDENTLY ADJUSTABLE HANDRAILS	•	•	○	N/A
PATIENT SAFETY BELT	•	•	•	○
EXERCISE BAND ATTACHMENT POINTS	•	•	•	N/A
NUMBER OF EXERCISE BAND ATTACHMENT POINTS	6	6	6	N/A
VARIABLE STABILITY PLATFORM	•	•	•	N/A
RESISTANCE MECHANISM	AIR PRESSURE (PSI)	AIR PRESSURE (PSI)	AIR PRESSURE (PSI)	FORCE PLATE
PLATFORM ADJUSTABILITY	21 LEVELS (0-20)	21 LEVELS (0-20)	21 LEVELS (0-20)	NONE
VERTICAL PLATFORM TILT RANGE	+/- 20°	+/- 20°	+/- 20°	NONE
HORIZONTAL PLATFORM TILT RANGE	360° ROM	360° ROM	360° ROM	NONE
COMPUTERIZED MECHANICAL PUMP	•	•	•	N/A
BALANCE BASELINE ASSESSMENT	•	•	•	•
COGNITIVE PARKINSON'S SOFTWARE	•	○	○	○
LIMITS OF STABILITY TEST (LOS)	•	•	•	•
mCTSIB TEST	•	•	•	•
SCORE RANGE	0 - 5000+	0 - 5000+	0 - 5000+	0 - 5000+
GRAPH REPORTING	•	•	•	•
ANTERIOR/POSTERIOR RATIOS	•	•	•	•
3D INTERACTIVE GAMES	•	•	•	•
NUMBER OF 3D GAMES	8	8	8	2
ABILITY TO BACKUP PATIENT RECORDS	•	•	•	•
PORTABLE UNIT	○	○	○	•
PATIENT WEIGHT CAPACITY	500 LBS. (227 KG)	500 LBS. (227 KG)	500 LBS. (227 KG)	330 LBS. (150 KG)
PRODUCT DIMENSIONS	44" W x 58" L x 81" H	44" W x 58" L x 81" H	44" W x 58" L x 81" H	27" W x 37" L x 67" H
PRODUCT WEIGHT	315 LBS (143 KG)	315 LBS (143 KG)	315 LBS (143 KG)	112 LBS (50 KG)
PLATFORM WEIGHT	N/A	N/A	N/A	31 LBS (14 KG)
ESTIMATED SHIPPING DIMENSIONS	63" X 45" X 81"	63" X 45" X 81"	63" X 45" X 81"	47" X 27" X 37"
ESTIMATED SHIPPING WEIGHT	420 LBS. (191 KG)	420 LBS. (191 KG)	385 LBS. (175 KG)	75 LBS. (34 KG)
CRATE OPTION FOR INTERNATIONAL SHIPMENTS	•	•	•	○
ESTIMATED CRATED DIMENSIONS	45" W x 79" H x 61" L	45" W x 79" H x 61" L	45" W x 79" H x 61" L	N/A
ESTIMATED CRATED WEIGHT	620 LBS (281 KG)	620 LBS (281 KG)	620 LBS (281 KG)	N/A
OPTIONAL SHIPPING CRATE COST	\$395.00 USD	\$395.00 USD	\$395.00 USD	N/A

- Feature Included
- Feature NOT Included
- N/A Not Applicable

In the interest of product development and improvement, specifications are subject to change without notice.

- Feature Included
- Feature NOT Included
- N/A Not Applicable

FEATURES:	KORESWAY MOBILE	KOREBALANCE ACTIVE	KORESWAY MINI	KOREBALANCE TRAINER
ITEM #	KB MOBILE	KB ACTIVE	KB MINI	KB TRAINER
COMPUTERIZED BALANCE ASSESSMENT	•	•	•	○
SCREEN SIZE	17.3	17.3	17.3	N/A
TOUCHSCREEN	•	•	○	N/A
MEDICAL GRADE TOUCHSCREEN	○	○	○	N/A
FULLY ADJUSTABLE SCREEN	○	○	○	N/A
SCREEN ADJUSTABILITY	○	○	○	N/A
PLATFORM SIZE	23" DIAMETER	23" DIAMETER	23" DIAMETER	23" DIAMETER
STEP UP HEIGHT	2 ¾ "	4 ½"	2 ¾ "	4 ½"
PATIENT HANDRAILS	○	○	○	N/A
VARIABLE STABILITY PLATFORM	N/A	•	N/A	•
RESISTANCE MECHANISM	FORCE PLATE	AIR PRESSURE (PSI)	FORCE PLATE	AIR PRESSURE (PSI)
PLATFORM ADJUSTABILITY	NONE	16 LEVELS (0-15)	NONE	21 LEVELS (0-20)
VERTICAL PLATFORM TILT RANGE	NONE	+/- 20°	NONE	+/- 20°
HORIZONTAL PLATFORM TILT RANGE	NONE	360° ROM	NONE	360° ROM
COMPUTERIZED MECHANICAL PUMP	N/A	○	N/A	○
MANUAL HAND PUMP	N/A	•	N/A	○
BALANCE BASELINE ASSESSMENT	•	•	•	○
COGNITIVE PARKINSON'S SOFTWARE	○	○	○	○
LIMITS OF STABILITY TEST (LOS)	•	•	•	N/A
mCTSIB TEST	•	•	•	N/A
SCORE RANGE	0 - 5000+	0 - 5000+	0 - 5000+	N/A
GRAPH REPORTING	•	•	•	N/A
ANTERIOR/POSTERIOR RATIOS	•	•	•	N/A
3D INTERACTIVE GAMES	•	•	•	N/A
NUMBER OF 3D GAMES	2	2	2	N/A
ABILITY TO BACKUP PATIENT RECORDS	•	•	•	N/A
PORTABLE UNIT	•	•	•	•
PATIENT WEIGHT CAPACITY	330 LBS. (150 KG)	500 LBS. (227 KG)	330 LBS. (150 KG)	500 LBS. (227 KG)
PRODUCT DIMENSIONS	23" W x 50" L x 65" H	23" W x 50" L x 65" H	23" W x 50" L x 65" H	23" W x 50" L x 65" H
PRODUCT WEIGHT	41 LBS (7 KG)	48 LBS (22 KB)	15 LBS (7 KG)	15 LBS (7 KG)
PLATFORM WEIGHT	15 LBS (14 KG)	23 LBS (10) KG)	9 LBS (14 KG)	22 LBS (10 KG)
ESTIMATED SHIPPING DIMENSIONS	27" X 27" X 9" 29" X 22" X 10" 24" X 7" X 5"	27" X 27" X 9" 28" X 27" X 9"	27" X 27" X 9"	27" X 27" X 9"
ESTIMATED SHIPPING WEIGHT	20 LBS (9 KG) 26 LBS (12 KG) 6 LBS (3 KG)	26 LBS (12 KG) 30 LBS (14 KG)	18 LBS (8 KG)	25 LBS. (11 KG)
CRATE OPTION FOR INTERNATIONAL SHIPMENTS	○	○	○	○